

Engineered for **performance**

Power. Precision. Productivity.



Mass spectrometry

TSQ Triple Quadrupole MS Portfolio

Quantify what is important for you, day after day, without compromise

Overview

The Thermo Scientific™ TSQ™ Triple Quadrupole Mass Spectrometer (MS) portfolio—Thermo Scientific™ TSQ Quantis Plus™ MS, Thermo Scientific™ TSQ Altis Plus™ MS, and Thermo Scientific™ TSQ Certis™ MS—provides a scalable range of triple quadrupole mass spectrometers designed to meet the needs of simple routine to highly challenging quantitative applications. Built from four decades of leadership in triple quadrupole mass spectrometry, each system delivers results.

From industrial, food safety, and environmental testing to clinical research and biopharma, the TSQ mass spectrometry portfolio delivers confident quantitation—no matter the application or the challenge.

Key benefits

- **Scalable performance** from routine to high-end systems
- **Consistent workflows and software** across all platforms
- **Advanced automation** for increased productivity

Essential specifications

	TSQ Quantis Plus MS	TSQ Altis Plus MS	TSQ Certis MS
			
Mass range (Da)	2–3,000	2–2,010	2–2,010
Scan speed (Da/s)	20,000	20,000 (H2 2026)	20,000
Acquisition speed (SRMs/s)	800	800	>900
Polarity switching (ms)	5	5	<5
Resolution (FWHM)	0.4	0.2	0.2
Linear dynamic range	>10 ⁶	>10 ⁶	>10 ⁶
Minimum SRM dwell (ms)	0.25	0.25	0.25
IDL Reserpine (HESI)	1.25 fg on column; IDL < 0.50 fg	0.5 fg on column; IDL < 0.20 fg	0.4 fg on column; IDL < 0.20 fg
IDL Chloramphenicol (HESI)	3.5 fg on column; IDL < 1.50 fg	0.6 fg on column; IDL < 0.25 fg	0.5 fg on column; IDL < 0.25 fg

Essential specifications			
	TSQ Quantis Plus MS	TSQ Altis Plus MS	TSQ Certis MS
HESI / APCI ion source (max °C)	Thermo Scientific™ OptaMax™ NG (550 °C)	Thermo Scientific™ OptaMax™ NG (550 °C)	Thermo Scientific™ OptaMax™ Plus (650 °C)
Compatible Thermo Scientific™ ion sources	• OptiSpray™ • EASY-Spray™ • Nanospray Flex™ • VeriSpray™ • OptaMax™ Duet • FAIMS Pro Duo	• OptiSpray™ • EASY-Spray™ • Nanospray Flex™ • VeriSpray™ • OptaMax™ Duet • FAIMS Pro Duo	• OptiSpray™ • EASY-Spray™ • Nanospray Flex™ • VeriSpray™ • FAIMS Pro Duo
Ion flight path	• Ion transfer tube • Stacked ring ion guide (SRIG) • Ion beam guide with neutral blocker • Thermo Scientific™ ARC™ (Active Reaction Collision) II collision cell with axial DC field • QR4 segmented quadrupole mass filter with hyperbolic surfaces	• High-capacity transfer tube • Electrodynamical ion funnel (EDIF) • Ion beam guide with neutral blocker • Thermo Scientific™ ARC™ (Active Reaction Collision) II collision cell with axial DC field • QR5 Plus segmented quadrupole mass filter with hyperbolic surfaces	• High-capacity transfer tube • Electrodynamical ion funnel (EDIF) • Q00 ion guide with mass filtering capability • Q0 ion beam guide with neutral blocker and axial DC field • Thermo Scientific™ ARC™ (Active Reaction Collision) II collision cell with axial DC field • QR5 Plus segmented quadrupole mass filter with hyperbolic surfaces
Dimensions (d, w, h)	33 × 30 × 27 in (840 × 760 × 680 mm)	33 × 30 × 27 in (840 × 760 × 680 mm)	33 × 26 × 25.6 in (833 × 703 × 660 mm)
Predictive maintenance	Yes	Yes (H2 2026)	Yes

 Learn more at thermofisher.com/tsq

General Laboratory Equipment – Not For Diagnostic Procedures. © 2026 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. Thermo Scientific™ general lab equipment, consumables, and software are not intended for *in vitro* diagnostic purposes in accordance with our product documentation, manuals, and labels. They are designated for General Laboratory Use Only. Not for use in diagnostic procedures. These products have not been tested or validated for such applications and their use for *in vitro* diagnostic purposes may result in inaccurate results and, more seriously, health and safety risks. They are designated (designed, intended) for research, academic, and industrial purposes. **PS004485-EN 0726**